

Alliance e2695 System Specifications

Solvent management system

Number of solvents	1~4
Flow rate range	0.01~10.0mL/min (使用時 0.05~5.0mL/min)
Maximum operating pressure	5000psi (0.01~3.0mL/min)

Sample management system

Number of sample vials	120vials (5carousels x 24vials)
Sample loops	100mL
Sample injection volume	5~80mL
Sample temperature control	4~40°C

Other

Column heater	20°C (room temperatute +5°C) ~60°C
Detector	UV/VIS (190~800nm, Photodiode array equipped)